


LINCOLN Diagnostics INC.

UNITEST[®] PC

PAIN CONTROL



The UniTest PC **TOUCH POST** activates nerve cells to inhibit pain.



UNITEST PC FITS INTO A 40-WELL TRAY.




ALLOW TOUCH POST TO CONTACT SKIN FOR ONE SECOND THEN PRESS UNITEST PC FIRMLY INTO SKIN.



4x magnification

ALL TEST SITES SHOULD SHOW THIS IMPRESSION ON THE SKIN WITH PROPER PRESSURE.



NEG
CNTL

PROVIDES EASY-TO-READ AND HIGHLY REPRODUCIBLE HISTAMINE RESULTS.

Lincoln Diagnostics, Inc. has been innovating the allergy skin testing industry since introducing Multi-Test[®] more than 40 years ago.

Now we're taking on pain. Our **UniTest PC** alleviates the pain of skin testing using the neurological gate-control theory, which involves the activation of nerve cells that block pain transmission to suppress pain. UniTest PC offers high sensitivity, high specificity, and low variability in results—with virtually no pain.

- Provides well-defined wheals with very little tissue trauma.
- Dipwell Tray and lid protect allergens from contamination.
- Dipwell Trays are stackable for convenient storage.
- Manufactured under ISO 13485 standards and OSHA compliant.
- Made in the USA.



Scan the QR code to access our UniTest PC Knowledge Base, which includes published data, videos, updated practice parameters, and additional documentation.

LINCOLN
Diagnostics INC.

The Leader in Allergy
Skin Testing Devices

(800) 537-1336 • 217-877-2531

www.lincolndiagnostics.com

